



TECO DIAGNOSTICS
1268 N. Lakeview Ave.
Anaheim, CA 92807
1-800-222-9880

BILIRUBIN TOTAL TC MATRIX-240/480

INTENDED USE

The test is applied for the quantitative determination of total bilirubin in human serum and plasma.

SUMMARY AND EXPLANATION OF THE TEST

Bilirubin, in bile, is derived primarily from the breakdown of hemoglobin when senescent red blood cells are phagocytized. Normally, about 6 to 6.5 g of hemoglobin in aged red blood cells are broken down daily in an adult to form about 220 mg of bilirubin; another 50 to 60 mg of bilirubin originates from other sources. The exact mechanism to form bilirubin is not well understood. But it is known that the heme group of hemoglobin converts to bilirubin in reticuloendothelial system, which binds to albumin in plasma and esterified to bilirubin digluconide (BDG), which is then secreted from liver as a waste product. It is present in serum in the free and conjugated forms. An increase in the formation or retention of bilirubin in the body results in increased levels of serum Bilirubin and jaundice. This hyperbilirubinemia is classified as either pre-hepatic, hepatic or post-hepatic depending on the principal cause of condition. Therefore, determination of the total bilirubin and its conjugated (direct) bilirubin is important for the differential diagnosis of hyperbilirubinemia.

TEST PRINCIPLE (*Colorimetric diazo method*)

Bilirubin reacts with diazotized Sulphanilic acid in an acidic environment to form the red azobilirubin. The intensity of the resulting color is measured by absorbance reading at a wavelength of 505 nm and is proportional to the total bilirubin concentration in the sample. As with the Jendrassik-Gr6f method, caffeine is present in the reagent as a reaction accelerator. Benzoate also has a similar effect. Surfactant as a solubilizing agent is also present in the reagent.

Total Bilirubin + diazo + H⁺ → azobilirubin (Blue color)

REAGENT COMPONENTS

Reagent 1

Sodium benzoate : ≤ 0.30 M
Sodium acetate : ≤ 0.50 M
Caffeine : ≤ 0.15 M

Surfactant

Reagent 2

Sulphanilic acid : ≤ 0.12 M
Hydrochloric acid : ≤ 0.22 M

Reagent preparation

No preparation is required.

REAGENT STORAGE AND STABILITY

Reagents are stable at 2-8°C until the expiration date stated on the label which is only for closed vials.

Once opened vials are stable for 30 days at 2-8°C in optimum conditions. On board stability is strongly related to auto analyzers' cooling specification and carry-over values.

Reagent stability and storage data have been verified by using Clinical and Laboratory Standards Institute (CLSI) EP25-A protocol.

Total Bilirubin activity stability in serum and plasma: 1 days at 20-25 °C, 7 days at 2-8°C and 6 months at -20°C

SPECIMEN COLLECTION AND HANDLING

Serum and plasma can be used and are collected according to the standard procedures. For plasma, sample collection tubes with Li heparin, Na heparin, K2-EDTA or K3-EDTA must be preferred. Contact with light must be

avoided. Non-lipemic samples must be used. Multiple samples freezing and thawing should be avoided

MATERIALS NEEDED BUT NOT SUPPLIED WITH REAGENT KIT

Arcal Auto Calibrator

At least two levels of control material.

CALIBRATION AND QUALITY CONTROL

Calibration: The assay requires the use of an Arcal Auto Calibrator.

Calibration stability is 30 days. Calibration stability depends on the application characteristics and cooling capacity of the autoanalyzer used.

Control: Commercially available control material with established values determined by this method can be used.

At least two level controls must be run once in every 24 hours. Each laboratory should determine its own quality control scheme and procedures. If quality control results are not within acceptable limits, calibration is required.

REFERENCE INTERVALS / MEDICAL DECISION LEVELS

Adults : ≤ 0.40 mg/dL
Newborn : Up to 24 hours : < 6.0 mg/dL
Up to 48 hours : < 10.0 mg/dL
3-5 days old : < 8.0 mg/dL
7 days old : < 10 mg/dL

PROCEDURE

Settings for TC-Matrix 240/480

Test Name:	D-BIL	R1:	147
Full Name:	Bilirubin Total	R2:	32
Pri. Wave:	505 nm	SAMPLE VOLUME:	5.7
Sec. Wave:	700 nm	Calibration Type:	2 point linear
Assay/ Point:	2-point end	K Value:	/
Start - End:	9 - 22	Point:	2
Decimal place:	2	Blank Type:	Reagent
Unit:	mg/dL	Point 0 (Blank) Con.:	0.0
Linearity Range:	0.1000 - 20.000	Point 1 (STD) Con.:	Calibrator / standard
Correlation Factor:	1.0000 - 0.0000		

PERFORMANCE CHARACTERISTICS

Measuring Interval

According to CLSI EP34-ED1:2018, "Measuring Interval" refers to the interval where the analyte concentration is measured with intended accuracy in terms of medical and laboratory requirements without dilution, concentrating or any kind of pre-treatment that is between the analyte's lower limit of quantitation (LLoQ) and upper limit of quantitation (ULoQ).

The determined analytic measuring interval for Bilirubin Total is 0.15 – 20 mg/dL.

Detection Capability

Limit of Detection (LoD): 0.1 mg/dL

Limit of Quantitation (LoQ): 0.15 mg/dL

Note: LoQ values are based on Coefficient of Variation Percentage (CV) ≤ 20%. LoD and LoQ values have been verified by using CLSI EP17-A2:2012 protocol.

Linearity

This method shows measurement linearity in the activities up to 20 mg/dL. Autoanalyzer's auto-dilution system can be used if the concentrations have higher values. See device manual for further information.

For the manual dilution procedure, dilute the sample 1:5 using 0.90% isotonic. After this process, multiply the result of the reworked sample by the dilution factor. Do not report the sample result after dilution if it is marked as lower than the linear lower limit. Rerun with a suitable dilution. Linearity Studies data have been verified by using CLSI EP06-A:2003 protocol.

Precision

Running system has been developed according to 20x2x2 "The Single Site" protocol. Repeatability and Within-Laboratory Precision/Within-Device values have been obtained according to the running results.

According to the protocol in use, 2 separate runs per day have been made for 20 days (no obligation for consecutive days). This protocol has been applied to each low and high samples separately and 80 results have been obtained for each one. Statistically, the results have been obtained using 2-factor Nested-ANOVA model.

Table 1. Total Bilirubin Repeatability (Within Run) results obtained from samples in two different concentrations

Mean Concentration	SD	CV%	n
0.68 mg/dL	0.01	2.25	80
5.18 mg/dL	0.02	0.55	80

Note: This working system has been named "Within-Run Precision" in the previous CLSI - EP05-A2 manual.

Table 2. Total Bilirubin Repeatability (Day to Day) results obtained from samples in two different concentrations

Mean Concentration	SD	CV%	n
0.68 mg/dL	0.03	3.74	80
5.18 mg/dL	0.15	2.90	80

Note: This working system has been named "Total Precision" in the previous CLSI - EP05-A2 manual.

Method Comparison

As a result of the statistical evaluation of the method comparison data: Passing-Bablok equation: $y = 1.179x + 0.184$ mg/dL
 $r = 0.96$

Interference

Endogenous interferant and analyte concentrations that have been used in the Total Bilirubin scanning tests has been determined according to "CLSI EP37-ED1:2018" and "CLSI EP07-ED3:2018" manuals.

The total acceptable error rate, which is going to be used to detect whether the observed differential value obtained from Total Bilirubin interference scanning test is appropriate, is determined as $\pm 10\%$.

In Total Bilirubin test results, no significant interaction has been observed in the determined endogenous interferant and analyte concentrations or between interferants and analyte.

Interferant Concentration	Total Bilirubin Target (mg/dL)	N*	Observed Recovery %
Indicane (Indoxyl sulfate) 0.08 mmol/L	1.06	3	100

* Total acceptable error rate determined as interference limit and repeatability (within run) pre-detected for the related method were used for the calculations of how many times the control and test samples prepared as a serum pool are going to

be run repetitively. In the calculations, the accepted error rate for type 1 (α error) was 5% and for type 2 (β error) was 10% (90% power).²⁷

WARNINGS AND PRECAUTIONS

IVD: For in Vitro Diagnostic use only. Do not use expired reagents.

Reagents with two different lot numbers should not be interchanged.

For professional use.

Follow Good Laboratory Practice (GLP) guidelines.

CAUTION: Human source samples are processed with this product. All human source samples must be treated as potentially infectious materials and must be handled in accordance with OSHA (Occupational Safety and Health Administration) standards

REFERENCES :

1. Malloy, H.T. and Evelyn, K.A. 1960, J.Biol. Chem. 119, 48.
2. Valdez, O.S.et.al.1971 J.Pediatr 79, 1015.
3. Tietz, N.W., "Specimen Collection and Processing; Sources of Biological Variation," Textbook of Clinical Chemistry, 2nd Edition, W.B. Saunders, Philadelphia, PA (1994).
4. Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1999:1136-7
5. Winsten S, Cehelyk B. A, rapid micro diazo technique for measuring total bilirubin. Clin Chim Acta 1969;25(3):441-6
6. Jendrassik L, Grof P. Vereinfachte photometrische Methoden zur Bestimmung des Blutbilirubins. Biochem Z 1938;297:81-9.
7. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline. CLSI Document EP25-A. Wayne, PA: CLSI; 2009.
8. Bhutani, V. K., Johnson, L., & Sivieri, E. M. (1999). Predictive ability of a pre-discharge hour-specific serum bilirubin for subsequent significant hyperbilirubinemia in healthy term and near-term newborns. Pediatrics, 103(1), 6-14. <https://doi.org/10.1542/peds.103.1.6>
9. Clinical and Laboratory Standards Institute (CLSI). Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline - Third Edition. CLSI Document EP28-A3c. Wayne, PA: CLSI; 2010.
10. Clinical and Laboratory Standards Institute (CLSI). Establishing and Verifying an Extended Measuring Interval Through Specimen Dilution and Spiking - 1st Edition. CLSI Document EP34. Wayne, PA: CLSI; 2018.
11. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline - Second Edition. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
12. Clinical and Laboratory Standards Institute (CLSI). Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach - 1st Edition. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
13. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline - Third Edition. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
14. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
15. Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.
16. Clinical and Laboratory Standards Institute (CLSI). Supplemental Tables for Interference Testing in Clinical Chemistry - First Edition. CLSI Document EP37. Wayne, PA: CLSI; 2018.
17. Clinical and Laboratory Standards Institute (CLSI). Interference Testing in Clinical Chemistry - Third Edition. CLSI Document EP07. Wayne, PA: CLSI; 2018.
18. CLIA proficiency testing criteria for acceptable analytical performance, as printed in the Federal Register July 11, 2022;87(131):41194-242.

B576-400TC2/TC4: 05/2024



TECO DIAGNOSTICS
1268 N. LAKEVIEW AVE.
ANAHEIM, CA 92807
U.S.A.