



## TECO DIAGNOSTICS

1268 N. Lakeview Ave.  
Anaheim, CA 92807  
1-800-222-9880

# COPPER CONTROL LEVEL I-II

### INTENDED USE

Copper Control level I-II are for quality control of the Copper Liquid assay.

### Materials Required But Not Provided:

1. Class A volumetric pipette for liquid transfer
2. Distilled or deionized water meeting the specifications equivalent to USP (United States Pharmacopeial Convention) purified water.

### CONTROL STABILITY

Temperature-Conditions	Stability
Unopened at 2-8°C	Expiry date on the vial
Opened and stored at 2-8°C	45 days

### PREPARATION OF CONTROL

Control is ready for use.

### INDICATIONS OF INSTABILITY OR DETERIORATION

Presence of extreme turbidity or microbial growth may indicate deterioration.

### PRECAUTION

Human source material. Treat as potentially infectious material. Each plasma donor used in the preparation of this product has been tested by an FDA-approved method and found negative for the presence of HIV 1/2 HBsAg, HCV, HIV-Ag antibodies. However, none of the known testing methods can offer complete assurance that the hepatitis B virus, Human Immunodeficiency Virus (HIV) or infectious agents are not present. All human-based products should be handled in accordance with Good Laboratory Practice (GLP) principles using appropriate precautions.

### WARNINGS

IVD: For in Vitro Diagnostic use only.

Do not use expired reagents.

Reagents with two different lot numbers should not be interchanged.

For professional use.

Follow Good Laboratory Practice (GLP) guidelines.

**CAUTION:** Human source samples are processed with this product. All human source samples must be treated as potentially infectious materials and must be handled in accordance with OSHA standards.

### Danger

H317: May cause allergic skin reaction.

### Precaution

P280: Use protective gloves / clothes / glasses / mask.

P264: Wash your hands properly after using.

P272: Contaminated work clothes should not be allowed to be used outside of the workplace.

### Intervention

P302+P352: Wash with plenty of water and soap if it contacts with skin.

P333+P313: Seek medical help if it irritates your skin or develops rash.

P362+P364: Remove contaminated clothes and wash properly before using.

### Disposal

P501: Dispose the vials and contents according to the local regulations.

### REFERENCES

1. Burtis CA, Ashwood ER, Bruns DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed. St. Louis, MO, Elsevier Saunders; 2006:2263.

2. Occupational Safety and Health Standards: bloodborne pathogens. (29 CFR Part 1910.1030). Fed. Register. July 1, 2001; 17:260-273.

3. Directive 2000/54/EC. Official Journal of the European Communities No. L262 from October 17, 2000.

4. Tietz, N. W. (1983) Clinical guide to Laboratory Tests, p384 W.B. Saunders Co., Philadelphia.

**C200-CI4&CI4: 03/2024**

Manufactured by:



TECO DIAGNOSTICS  
1268 N. LAKEVIEW AVE.  
ANAHEIM, CA 92807  
U.S.A.